

# Institut La Claire EM2

## ACTIVE DRY YEAST

SACCHAROMYCES CEREVISIAE STRAIN ISOLATED AND SELECTED BY THE MASSEY UNIVERSITY IN NEW ZEALAND.

Conforms with the International Oenological Codex. Not derived from genetically modified organisms.  
Allergen free.

### Oenological properties

Distinguishing feature: high polysaccharide production. EM2 is a varietal strain, meaning that it is equipped with enzymes that are able to free and enhance the terpene aromatic components peculiar to each grape variety. It also boasts reliable fermentation marked by a high percentage of free SO<sub>2</sub>, as well as low acetaldehyde and volatile acidity levels. Thanks to the release of polysaccharides, it also helps boost the roundness and softness of well-developed wines.

### Advanced properties

The aromatic profile of EM2 is characterized by a positive and intriguing note reminiscent of the flatness of Sauvignon. It also has the right characteristics for producing high-range white wines, including those fined or fermented in barriques – their varietal characteristics are enhanced thanks to the fruity and floral notes enhanced by this strain. By releasing polysaccharides, this strain can also contribute to colour stabilization in red wines lending them greater intensity and stronger nuances. EM2 responds best when nourished with substantial quantities of RAN (readily assimilable nitrogen).

### Varietal white wines and beyond

Given its unique characteristics, EM2 is naturally suited to the production of strong, well-structured white wines such as Sauvignon, Chardonnay, Vermentino and Pinot Grigio; it also enhances Fiano, Müller-Thurgau and Malvasia. The roundness, softness and colour stability deriving from the polysaccharides also make this strain an interesting option for red wines such as Valpolicella, Teroldego and Cabernet Franc.

### Composition

Yeast, E491.

### Characteristics

Appearance: small rods.

Colour: light ochre.

Alcohol production: 15% v/v

Classification: *Saccharomyces cerevisiae*

Cell count: > 10bn. live cells/gram

Optimum temperature: 15-30 °C

### Dosage

10-25 g/hl for fermentation.

30-50 g/hl for stuck fermentation or in the most difficult cases.

### How to use

Rehydrate the yeast in a suitable, sanitized containers with at least 10 parts water at a temperature of around 40°C (from 30°C to 45°C). After rehydration, leave the mixture to settle for 10 minutes. Start stirring again, continue hydration for a further 10 minutes and then add to the must to be fermented straight away. For top performance, it is advisable to add the same proportion of the Ecobiol Pied de Cuve nutrient.

When using for the second fermentation, acclimatize the yeasts before inoculation.

### Storage

Store in a cool, dry environment.

Use by the date printed on the package.

Once open, store in a refrigerator at +4°C.

### Warnings

Do not rehydrate in must or in cold water. It is advisable not to leave the yeast in the water for any longer than the recommended time (do not exceed 30 minutes total).

Once the pack has been opened, it is preferable to use the whole content immediately.

### Pack sizes

Code 107291 – 500g vacuum packs



**Perdomini-IOC S.p.A.**

Via Salvo D'Acquisto, 2 - 37036 S. Martino Buon Albergo (VR) Italy

tel. +39-045-8788611 r.a. fax +39-045-8780322

fax uff. vendite +39-045-8780122

www.perdomini-ioc.com - info@perdomini-ioc.com



At the time of printing, the information contained in this publication has been checked and deemed correct.

Perdomini-IOC S.p.A. reserves the right to change the characteristics of its products without forewarning.